

USE MORE CELLS

For a free 6-pack of one of our many Thermo Scientific HyClone Classical Media products, visit the website below. Media products available with this offer listed below. No cash value. Limit one free 6-pack per customer.

www.thermoscientific.com/usemorecells

Thermo Scientific HyClone is dedicated to meeting your cell culture needs with the highest quality sera, media and reagents.

MEDIA TYPE	MODIFICATIONS	PART #
Dulbecco's Modified Eagle's Medium (DMEM)	with High Glucose, with 4.0 mM L-glutamine, without Sodium Pyruvate	SH30022.FS
	with High Glucose, without L-glutamine and Sodium Pyruvate	SH30081.FS
	with High Glucose, with 4.0 mM L-glutamine with Sodium Pyruvate	SH30243.FS
	with High Glucose, without L-glutamine, with Sodium Pyruvate	SH30285.FS
	with Low Glucose, 4.0 mM L-glutamine, and 110 mg/Sodium Pyruvate	SH30021.FS
DMEM/F12 1:1	with 2.5 mM L-glutamine, with HEPES	SH30023.FS
	without L-glutamine, with HEPES	SH30126.FS
	with L-glutamine, without HEPES	SH30271.FS
Ham's Nutrient Mixture	with 1.0 mM L-glutamine	SH30026.FS
Iscove's Modified Dulbecco's Medium (IMDM)	with L-glutamine, with HEPES, without Alpha-Thioglycerol	SH30228.FS
McCoy's Medium	McCoy's 5A with L-glutamine	SH30200.FS
Minimum Essential Medium (MEM)	with Earle's Balanced Salt Solution (EBSS), with 2.0 mM L-glutamine	SH30024.FS
	with Earle's Balanced Salt Solution (EBSS), without L-glutamine	SH30244.FS
	Alpha Modification, with L-glutamine, Ribo- and Deoxyribonucleosides	SH30265.FS
RPMI 1640 Medium	with 2.05 mM L-glutamine	SH30027.FS
	without L-glutamine	SH30096.FS
	with 25 mM HEPES, with L-glutamine	SH30255.FS

For a complete listing of all the Thermo Scientific HyClone quality products for cell culture, Visit www.thermoscientific.com/hyclone

To see the details and request a free 6-pack of one of the medias above, Visit www.thermoscientific.com/usemorecells